

JAMIA HAMDARD Department of Paramedical Sciences School Of Nursing Sciences & Allied Health Medical Laboratory Techniques

S. No.	, Equipment	Specification	Qty.
1.	Pipettes (2ml)	Volume -2ml	200
2.	Pipettes (5 ml)	Material - Silica gel Volume- 5ml	200
	Pipettes (10 ml)	Material - Silica gel Volume- 10 ml	25
3.		Material - Silica gel	
4.	Micropipettes fixed (10 ul)	Volume- 10ul Material -Polypropylene	
5	Micropipettes fixed (50 ul)	Volume - 50ul Material -Polypropylene	100
6	Micropipettes fixed (100 ul)	Volume -100ul Material -Polypropylene	
7	Micropipettes fixed (1000 ul)	Volume - 1000ul Material -Polypropylene	
8. Glass flask 250 ml		Volume - 250ml	20
		conical, Narrow Mouth Erlenmeyer, Graduated,	
9	Glass flask 500 ml	Volume - 500ml Material - Borosilicate glass	20
10	Glass flask 1000ml	Volume - 1000ml Material - Borosilicate glass	20
11	Beaker 250 ml	250ml Heat Capacity(Degree celsius) 560 Degree Celsius Material : Glass	20
12	Incubator	Power Consumption - 230 volts Dimensions - 18*18*18 Model no 8	
13.	Centrifuge	Speed - 4000RPM Motor- Brush less	20
14 Beaker 500 ml		LIDI C	

		capacity: 500ml	
15	Beaker 1000 ml	Capacity: 1000 ml Material: Borosilicate Glass	20
16	Colorimeter Lab Junction Colorimeter (Digital Photo Colorimeter), Model: LJ-312	Ean: 4379746680781 Item Weight: 1.00 kilograms Model Number: LJ-312	01
17	Microscope LED lamp		02
18	Refrigerator	Blue Glass Door Single Zone Temperature Model no - MTGDR11W Capacity - 500L	1
19	Spirit Lamp aluminium	Color: Silver Material: Aluminium Packaging Type: Plastic Wrapping Shape: Burner Transparent: No	25
20	RBC pipette	Brand:STANDARD Color: TRANSPERANT Material: GLASS	50
21	WBC pipette	Capacity: 1 ml Usage/Application: Used for mixing diluting fluid and blood Material: Glass Color: Transparent	50
22	Urinometer	For specific gravity of urine Analogue display Manually operated	10
23	Slide Racks	Polylab 18 Places PP 3-Tier Test Tube Stand for 25mm Tube Material: Polypropylene Item Code: 77710	20
24	Specimen slide for teaching	SKU: LAB. MIC.23068266 Type of Product: Microslides Plain Gr. Edges Model No: 9100P02 Sub Type: Double Frosted Ground Edges Dimension (mm): 76X26X1	20

25	Tube rack 13mm x 31 hole	Brand Name: Clear & Sure Item Weight: 200 grams Material: plastic Rack D18,40 Holes Size 18 MM	20
26	Pasteur Pipette (500/pack)	Capacity 3.00 milliliters Marking Colour Transparent Material polypropylene Part Number KTN400	25
27	Modified neubauer's chamber	thick crystal slide with the size of a glass slide. (30 x 70 mm and 4 mm thickness)	25
28	Sahli's Hemoglobin meter	Brand :AQRA Model AQRA Sahli's Haemometer,Prismatic (Square Tube),Planoparallel (Square Tube / Round Tube) Application: Clinical Measuring range: 0-20 gm/dl Type of Tube: Square	25
29	Tourniquet with buckle	Brand: Neutral Length: 16 Inch Color: Blue Pack Size:200 Piece Per Box	10
30	'Test tube racks	Color: Blue Material: Plastic Brand: Satya moulds Size: 50ML	25
31	Westerngren ESR Tube and stands	Brand:Thermo Fisher Material: Plastic Shape: Round Minimum Order Quantity: 50 Piece(s)	25
32	Wintrobe ESR tube and stands	Color: Transparent Capacity:200mm & 300 mm Pack Contain: 50Tubes	50

.

		Material:Glass	
3	Rubber bulb (medium size)	Material:Rubber Color: Red Shape: Round	50
4	Rubbber bulb (small size)	Material: Natural Rubber Brand: Rubbertron Capacity: 25 ML Diameter: 53 MM	50
5	Filter paper sheets	material:cellulose sheets packaging: pack of 100 manufacturer/tradename:Wha tman 1001-917 Whatman Article No. 28413937 parameter:0.25 psi wet burst 150 sec/100 mL speed (Herzberg) W × L:460 mm × 570 mm thickness: 180 µm pore size:11 µm (Particle retention) basis weight:87 g/m2	50
36	Test Tubes (5 ml)	Capacity: 5 ml Material: Plastic Color: Transparent	250
37	Test Tubes (12 ml)	Brand: Recombigen Volume: 12 ml Packaging Size: 100 pcs Cap:Screw Cap Material: Polypropylene Color:Transparent	250
38	Water bath	Display Mode: Analog Level of Automation: Semi Automatic Frequency: 33 KHz Brand: Samarth Electronics Material: Stainless Steel Digital Timer: 20 Min	01
39	9 Glass slide	Size: 75mm x 25mm Material: Glass Thickness Of Slide: 1.25mm	500

10	Chromatography chamber	Brand: Labline Material: Glass, Wood Shape: Cuboid	01
41	Leishman staining kit	Packaging Size: 25 kg Usage/Application: Leshmans Stain BIOLOGICAL STAINS	20
42	Antihuman globulin AHG	5ml Biotechnika	10
43	ABO Reagent	Product Volume-5ml Storage condition-2-8 degree celcius	10
44	AFB Staining kit	Product Volume-50gm Storage condition -room tempertaure	10
45	Sodium nitroprusside	Product Volume-500ml Storage condition-2-8 degree celcius	10
46	Ammonium sulphate	Product Volume-500ml Storage condition-2-8 degree celcius	10
47	Acetone	Product Volume-500ml Storage condition-2-8 degree celcius	10
48	Sulphur powder	Product Volume-150gm Storage condition -room tempertaure	10
49	Normal saline	Product Volume-50ml Storage condition-2-8 degree celcius	10
50	Gram staining kit	Product Volume-100gm Storage condition -room tempertaure	10
51	Albert staining kit	Product Volume-100gm Storage condition -room tempertaure	10
52	Gloves	Non powder , Latex Size-7.5	· 250 box
53	Blood collection vial (Red, green, voilet)	3ml each vial	250 box
54	Cotton	Type: Plain Sterility: Non-Sterile	10
55	Tube holder with vacutainer	Plastic Stand of Stx	10
56	Tissue paper roll	For clean slides and tubes before any investigation starts.	2

57	Sodium hydroxide	Product Volume-100gm Storage condition -room	2
		tempertaure	
58	Sodium chloride	Product Volume-100gm	2
		Storage condition -room	2
		tempertaure	
59	Sucrose '	Product Volume-100gm	2
		Storage condition -room	2
		tempertaure	
60	pH meter strip	Portable meters accucartely	2
		measure pH,conductivity,ion	-
		concentration, ORP	
61	Laminar air flow	Vertical laminar flow to	20
		protect sample, free flow of	20
		air.	
62	Drabkin's Reagent	Drabkin's reagent 100ml is	01
		used in measuring	01
		haemoglobin from blood	
		samples.	
63	Pipette stand	A pipette stand is a holder	25
	A STATE OF THE PARTY OF THE PAR	designed to provide safe and	23
		convenient storage for	
		manual and electronic	
		pipettes	
64	Coverslip	cover slip is to keep solid	300
		specimens pressed flat, and	
		liquid samples shaped into a	
		flat layer of even thickness	
65	M	(5mm)	
65	Microtrips (100ul)	Microtips are thin, high	500
		intensity probes which are	
		designed for processing small	
66	Mignatula (1000)	sample volumes	
00	Microtrips (1000ul)	Microtips are thin, high	500
		intensity probes which are	
		designed for processing small	
67	Cnatula	sample volumes	
07	Spatula	spatulas are used to transfer	25
		samples or chemicals from	
		their original containers to	
		weighing paper, weighing	
		boats, weighing bottles,	
		weighing funnels, or other	
		vessels or containers for	
58	Reagent bottle 500 ml	weighing.	
	heagent bottle 500 ml	Reagent bottles are	20
		excellent for storing	
		powders and liquids.	
		Reagent bottles in tinted	
		amber or red protect light-	

		111/	
		sensitive contents from UV	
		light, visible light and	
		infrared radiation.	
69	Reagent bottle 100ml	Reagent bottles are	20
09	(Amber colour)	excellent for storing	
	(Amber colour)	powders and liquids.	
		Reagent bottles in tinted	
	,	amber or red protect light-	
		sensitive contents from UV	
		light, visible light and	
		infrared radiation.	
		A glass stirring rod, glass rod,	10
70	Glass rod	stirring rod or stir rod is a	
		piece of laboratory equipment	
		used to mix chemical	
			25
71	Coplin jar	-coplin jar molded in PP can	23
		hold 10 slides of 1" x 3" size.	
		This jar is grooved inside in	
		such a way that it can hold	
		the slides erect.	25
72	Graduating measuring	A Measuring Cylinders is a	25
	cylinder (100 ml)	piece of laboratory equipment	
		that is use to measure the	
		volume of a any liquids,	
		chemicals or solutions during	
		the laboratory	
73	Hot plates	in laboratory settings, hot	02
13	Hot places	plates are generally used to	
		heat glassware or its contents.	
74	Forceps	primary purpose is to grasp,	50
/4	Гогсерз	retract, or stabilize tissue.	
75	Cell counter	Automated machine to	10
13	Cell codifice	calculate the Complete blood	
		count	
		Count	
76	Petri dish	A Petri dish is a shallow	25
10	retraisir	transparent lidded dish that	
		biologists use to hold growth	
		medium in which cells can be	
		cultured	
	C : Warratubo	The main aim of the capillary	500
77	Capillary tube	tube is to reduce pressure. It	
		is used in laboratories for	
		collecting the different	
		samples of liquid.	
		An Urinometer is a simple	25
78	Urinometer	piece of equipment for	
		determining urine specific	
		gravity	

79	Dipstick	A urinalysis is a test of your	50
		urine. It's used to detect and	50
		manage a wide range of	
		disorders, such as urinary	
		tract infections, kidney	
80	C	disease and diabetes.	
80	Sensitivity disk	These disks are used to test	25
	,	the effectiveness of particular	25
		antimicrobial agents	
		(antibiotics) against specific	
81	Platelet diluting and	target organisms.	
82	Platelet diluting fluid	used in hematology.	10
02	Syringes	A syringe is a small tube with	
		a thin hollow needle at the	50 box
83	Needlo (green 11)	end. (5ml)	
84	Needle (green and blue)	(green and blue)23g	50 box
04	Semi-Auto analyser	Semi-automated analysers are	
		faster, more precise chean	01
		and compact as compared to	
		Tully automated analysers and	
		are capable to perform tests	
		on whole blood, serum	
		plasma, cerebrospinal fluid	
35	Haematology A1	and urine samples	
	Haematology Analyser	Flaematology analysers are	0.1
		used to run tests on blood	
-		samples.	



Jamia Hamdard Department of Paramedical Sciences School of Nursing Sciences and Allied Health

List of Equipment's For BCLT

Equipment	Quantity	Specifications	
Blood Pressure Machine(Sphygmomanom eter)	3	Conventional Mercurial type BP Instrument	
2. 12 Channel ECG Machine	1	12 Lead ECG Machine with ECG Graph papers	
3. Multi Brand Echo Machine Portable	1	Multi Brand/Machine condition refurbished/Image Mode 2 D/Doppler Mode: Continuous Wave, Pulse Wave Doppler, Color Doppler Velocity, Power Doppler	
4. Ambu bag	5	Medical Rebreathing Bag/Ambu bag for Adult and Kids/Silicon Reusable Oxygen Reservoir Bag for Hospital, Clinic and Home(2600 ml)	
5. CPR Dummy Single Ultra- light CPR Training Manikin	1	Full Body CPR Dummy with monitor/portable plastic bag /CPR Operation Pad/CPR Face shield 50 pieces per Box/Lung Bag/Exchangeable skin /Guide Manual	
6. Stethoscope	5	Stethoscope for clinical examination/outer diameter 10mm/Inner diameter 4.8 mm/Tube Length 560	

		mm/sensitivity 3.2 to 26 DB from 50Hz to 100 Hz for Cardiology
7. 3D Heart model	1	Prosthetic Human heart Anatomical model, 2 parts, Showing 4 chambers, Valves, Major Blood vessels with detailed anatomical structures
8. Human Skeleton	2	Prosthetic Human Anatomical Skeleton, Medical Teaching Model/White



Jamia Hamdard Department of Paramedical Sciences School of Nursing Sciences and Allied Health New Delhi-110062

Equipments for the Optometry Program

S.N	Name of Equipment	Quantity	Specifications
1.	Spectacle lens grinding & surfacing machine semiautomatic for spherical lens and cylindrical lenses & Tools for grinding optical lenses of various curvatures gauge for measuring curvature of optical lenses.	01	Semiautomatic/Manual spherical & cylindrical spectacle lens grinding Machine, for giving power to glass lenses. One table for spherical lens power with half HP motor, one table for cylinder lens power with one HP motor & one table for polishing with aluminium Tub one HP power, Also accessories in the form of grinding tools of various curvatures.
2.	Manual lens edging machine with wheels	01	Manual hand edger with diamond wheels medium size, automatic water circulation, motorized
3.	Optical lens fitting set of callipers, screw drivers & Lens cutter	01 each/01 set	For adjustment of frames & lenses in a rimless, supra or full fitting stainless steel frame.
4.	Lens marking pencil	10	Black or blue colour for marking on glass
5.	Trial lens Box Set (Refraction set)	03	Lenses with Plano back vertex fitted in aluminium mounds, plus & minus lenses with accessories Maddox rod, occlude, blank, prisms, red and green glass, pin hole and stenopaic slit, cross cylinder +/- 0.25D +/- 0.50D and or PD ruler.
6.	Trial Frames	10	Metal adjustable,3 grooves, adjustable PD, regular model, light weight
7.	Vision Drum(Distance)	05	Snellen's chart: motorized with remote? Manual, with duochrome test, operable at 6 mts and or 3 mtrs ETDRS chart (Lovie-Bailey chart), Sherdan-garnier chart for children
8.	Near Vision Charts	05	Illuminated, plastic screens with metal frames for testing near vision.
9.	Axis marker	02	For marking asis on the finished spectacle lens
10.	Streak Retinoscope	04	Hiene's Model
11.	Glass blanks (unifocal) (consumable)	100	For grinding power to form finished spectacle lens
12.	Glass blanks (Bifocal)(consumable)	100	For grinding power to form finished spectacle lens



Jamia Hamdard Department of Paramedical Sciences School of Nursing Sciences and Allied Health , New Delhi-110062

Equipments for the Optometry Program

S.N	Name of Equipment	Quantity	Specifications
1.	Spectacle lens grinding & surfacing machine semiautomatic for spherical lens and cylindrical lenses & Tools for grinding optical lenses of various curvatures gauge for measuring curvature of optical lenses.	01	Semiautomatic/Manual spherical & cylindrical spectacle lens grinding Machine, for giving power to glass lenses. One table for spherical lens power with half HP motor, one table for cylinder lens power with one HP motor & one table for polishing with aluminium Tub one HP power, Also accessories in the form of grinding tools of various curvatures.
2.	Manual lens edging machine with wheels	01	Manual hand edger with diamond wheels medium size, automatic water circulation ,motorized
3.	Optical lens fitting set of callipers, screw drivers & Lens cutter	01 each/01 set	For adjustment of frames & lenses in a rimless, supra or full fitting stainless steel frame.
4.	Lens marking pencil	10	Black or blue colour for marking on glass
5.	Trial lens Box Set (Refraction set)	03	Lenses with Plano back vertex fitted in aluminium mounds, plus & minus lenses with accessories Maddox rod, occlude, blank, prisms, red and green glass, pin hole and stenopaic slit, cross cylinder +/- 0.25D +/- 0.50D and or PD ruler.
6.	Trial Frames	10	Metal adjustable,3 grooves, adjustable PD, regular model, light weight
7.	Vision Drum(Distance)	05	Snellen's chart: motorized with remote? Manual, with duochrome test, operable at 6 mts and or 3 mtrs ETDRS chart (Lovie-Bailey chart), Sherdan-garnier chart for children
8.	Near Vision Charts	05	Illuminated, plastic screens with metal frames for testing near vision.
9.	Axis marker	02	For marking asis on the finished spectacle lens
10.	Streak Retinoscope	04	Hiene's Model
11.	Glass blanks (unifocal) (consumable)	100	For grinding power to form finished spectacle lens
12.	Glass blanks (Bifocal)(consumable)	100	For grinding power to form finished spectacle lens

13.	Emeries powder of different size No.	4 packet each	For grinding power to form finished spectacle lens
	1.Carborandom grain		
	size 0.3mm.		
	2.Aluminium oxide		
	202 & 203 grain size ,		
	0.075mm & 0.01mm.		
	3.Iron Oxide grain		
	size 0.002mm		
1.4	(consumable)	10	For grinding power to form finished spectacle lens
14.	Cloth for grinding	10 meters	For grinding power to form missied spectacle lens
	(consumable)/polishin g pad		
15.	Cross cylinder 0.25-	02	Jackson's cross cylinder with handle & powers +/- 0.25D & +/-
15.	0.5D powers	02	0.5D
16.	Colour vision Ishihara	02	With 39 plates original Japanese for testing colour vision
10.	chart	02	With 37 places original supuness for costing second
17.	Hertel'sexophthalmo	01	For measuring exophthalmose
17.	meter	O1	To measuring exceptional
18.	PD Ruler	02	For measuring the pupillary distance
19.	Titmus Fly test	01	Stereopsis measuring for Adults & children with the symbols,
1),	Titings Lij test		stereo fly with stereopsis measuring,+/-4800 to 20 seconds of
			arc.
20.	Maddox wing	01	For measuring heterophoria, light in weight metal
21.	Diploscope with	01	For measuring diplopia, metal frame
	picture card		
22.	Placido' disc	01	Plastic hand held for seen corneal surface
23.	Eye model	05	For refraction training
24.	Frame heater	01	For heating of non-metal spectacle frames to fit lenses
25.	Radian	01	Measuring finish lens curvature -0 to 20 mm .reading -0.25
	Apparatus/Geneva lens measure		measure, scope +-20
26.	Low Vision kit	01	To teach students about low vision
27.	Geneva Lens Measure	01	For measuring the base curve of a lens
28.	Soft lenses, Toric	02	For contact lens fitting
	lenses, and RGP lens		
	trial set		
29.	Automated lensometer	01	For measuring the power of glasses
	(grand Seiko)		il and hoom onlitter
30.	Slit lamp with moving	01	Zeiss type slit lamp three steps with camera and beam splitter
	table		
21		0.1	Appasamy Synoptophore ASP 902, 430 mm (H) X 550 mm (c
31.	Synaptophore	01	265 mm (B).
20	Dist HOTY falling	01	Pediatric eye chart featuring HOTV optotypes for
32.	Dist.HOTV folding	01	screening/testing the vision of children and adults 3 years and
			older. This eye chart folds in half for added mobility or easy
			storage. Front side of eye chart features HOTV characters in
			proportionally-spaced lines (ETDRS format) for monocular and
			binocular (left eye / right eye) vision testing.
	Accommodation		Test and train for accommodative facility, accommodative rock,
33	i / iccomminuodation		
33.	Flippers Set	05	positive/negative relative accommodation and convergence and

	Trans a Ci		
34.	ETDRs Chart	01	For testing at 4 meters (13 ft.) Acuity range 40M to 2M (20/200, 6/60 to 20/10, 6/3 equiv.) Only standardized format accepted in studies sponsored by the National Eye Institute Designed for use in clinical studies and low vision evaluations where accurate, repeatable visual acuity measurements are required Lines are of equal difficulty Geometric progression in letter size from line to line Letter size is only variable Measures 62 X 65 cm (24.5 X 25.5 in.) Chart fits the ETDRS Illuminator Cabinet CAT. NO. 2425 Can be wall-mounted
35.	Pelli-Robson Contrast Sensitivity Chart	01	Features Traditional Sloan Optotypes For use at 100cm (1m) Includes two charts which may be wall or easel mounted and score sheets. Illumination of 85 cd/m² is key to good results Designed by DG Pelli, JG Robson and A J Wilkins, because they see contrast sensitivity as a realistic assessment of how well the patient sees large faint objects. The chart follows the luminance, font, and letter spacing recommendations of the Committee on Vision of the National Academy of Sciences and National Research Council (Adv Ophthal 41, 103-148, 1980). The letter sequences are organized into groups of three (triplets) with two triplets per line. Within each triplet, all letters have the same contrast. The contrast decreases from one triplet to the next, even within one line. The testing of contrast sensitivity on a regular basis is key to detecting early signs of diseases such as Parkinson's. Glaucoma, and Cataracts. By testing with several different contrast levels configured together, the patient's contrast sensitivity level is better assessed and records will be more precise.

1.3 Onkumun

... appendent itils

Jamia Hamdard Department of Paramedical Sciences School of Nursing Sciences & Allied health

B.Sc. Trauma Care Management

Equipment Required for Trauma Care Management

Endotracheal Tube 4 trachea for emergency ventilation Sizes-6,6.5,7,7.5,8 With curved blade (macintosh), straight blade (millers), Device used to examine the larynx for intubation Cylindrical Flexible rubber instrument used for making passages Board used for to make surface firm for effective CPR Used for lifting and moving the patient in case of trauma with starps to tie the patient(orange in color with grooves cut) Kit used for spine immobilisation Kit Splints (Rigid and Soft) AMBU Bag CPR Mannequin Device used for immobilizing broken bone of the body to prevent further trauma Device used for artificial ventilation (Sizes-adult, paediatrics) to learn how to give cpr (Brand- Laedral) Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Cannula used to place in the central veins Used to open a port in blood vessel for medication and sampl withdrawal (Sizes-16,8,20,22,24)	Equipment's	Quantity	Specification
Tube 4 Sizes-6,6.5,7,7.5,8 With curved blade (macintosh), straight blade (millers), Device used to examine the larynx for intubation Cylindrical Flexible rubber instrument used for making passages Bougie 2 Board used for to make surface firm for effective CPR Used for lifting and moving the patient in case of trauma with starps to tie the patient (orange in color with grooves cut) Spine Immobilisation Kit 2 Splints (Rigid and Soft) 4 Device used for immobilizing broken bone of the body to prevent further trauma Device used for artificial ventilation (Sizes-adult, paediatrics) To be fibrillator 1 Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Central Line 2 Cannula used to place in the central veins Used to open a port in blood vessel for medication and sample withdrawal (Sizes-16,18,20,22,24)			Tube which is passed through
With curved blade (macintosh), straight blade (millers), Device used to examine the larynx for intubation Cylindrical Flexible rubber instrument used for making passages Board used for to make surface firm for effective CPR Used for lifting and moving the patient in case of trauma with starps to tie the patient(orange in color with grooves cut) Spine Immobilisation Kit Splints (Rigid and Soft) AMBU Bag CPR Mannequin Defibrillator Central Line Cannula (diff size) Subject to sea to examine the larynx for intubation Cylindrical Flexible rubber instrument used for making passages Board used for to make surface firm for effective CPR Used for lifting and moving the patient in case of trauma with starps to tie the patient(orange in color with grooves cut) Kit used for spine immobilization with head blocks Ait used for immobilizing broken bone of the body to prevent further trauma Device used for artificial ventilation (Sizes-adult, paediatrics) to learn how to give cpr (Brand- Laedral) Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Cannula used to place in the central veins Used to open a port in blood vessel for medication and sampl withdrawal (Sizes-16,18,20,22,24)	Endotracheal		
Laryngoscopes 2	Tube	4	Sizes-6,6.5,7,7.5,8
Laryngoscopes 2 (millers), Device used to examine the larynx for intubation Cylindrical Flexible rubber instrument used for making passages Board used for to make surface firm for effective CPR Used for lifting and moving the patient in case of trauma with starps to tie the patient (orange in color with grooves cut) Spine Immobilisation Kit Splints (Rigid and Soft) AMBU Bag CPR Mannequin Defibrillator Central Line Cannula (diff size) Spine Immobilisation Cannula (diff size) Cannula (diff size) (millers), Device used to examine the larynx for intubation Cylindrical Flexible rubber instrument used for making passages Board used for to make surface firm for effective CPR Used for lifting and moving the patient in case of trauma with starps to tie the patient (orange in color with grooves cut) Kit used for spine immobilization with head blocks A bevice used for immobilizing broken bone of the body to prevent further trauma Device used for artificial ventilation (Sizes-adult, paediatrics) to learn how to give cpr (Brand-Laedral) Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Cannula used to place in the central veins Used to open a port in blood vessel for medication and sampl withdrawal (Sizes-16,18,20,22,24)			With curved blade
Laryngoscopes 2 the larynx for intubation Cylindrical Flexible rubber instrument used for making passages Board used for to make surface firm for effective CPR Used for lifting and moving the patient in case of trauma with starps to tie the patient(orange in color with grooves cut) Spine Immobilisation Kit 2 Device used for immobilizing broken bone of the body to prevent further trauma Device used for artificial ventilation (Sizes-adult, paediatrics) to learn how to give cpr (Brand-Laedral) Defibrillator 1 Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Cannula used to place in the central veins Used to open a port in blood vessel for medication and sample withdrawal (Sizes-16,18,20,22,24)			(macintosh), straight blade
Cylindrical Flexible rubber instrument used for making passages Board used for to make surface firm for effective CPR Used for lifting and moving the patient in case of trauma with starps) Spine Immobilisation Kit Splints (Rigid and Soft) AMBU Bag CPR Mannequin Defibrillator Central Line Cannula (diff size) CPR Board CPR Board CPR Board CPR Board I Cylindrical Flexible rubber instrument used for making passages Board used for to make surface firm for effective CPR Used for lifting and moving the patient (orange in color with grooves cut) Kit used for spine immobilization with head blocks Device used for immobilizing broken bone of the body to prevent further trauma Device used for artificial ventilation (Sizes-adult, paediatrics) to learn how to give cpr (Brand-Laedral) Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Cannula used to place in the central veins Used to open a port in blood vessel for medication and sample withdrawal (Sizes-16,18,20,22,24)			(millers), Device used to examine
Bougie 2 passages CPR Board 1 Board used for to make surface firm for effective CPR Used for lifting and moving the patient in case of trauma with starps to tie the patient (orange in color with grooves cut) Spine Immobilisation Kit 2 Splints (Rigid and Soft) 4 Device used for immobilizing broken bone of the body to prevent further trauma Device used for artificial ventilation (Sizes-adult, paediatrics) To learn how to give cpr (Brand-Laedral) Defibrillator 1 Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Central Line 2 Cannula used to place in the central veins Used to open a port in blood vessel for medication and sample withdrawal (Sizes-16,18,20,22,24)	Larvngoscopes	2	the larynx for intubation
Bougie 2 passages CPR Board 1 firm for effective CPR Used for lifting and moving the patient in case of trauma with starps to tie the patient(orange in color with grooves cut) Spine Immobilisation Kit 2 Splints (Rigid and Soft) 4 Device used for immobilizing broken bone of the body to prevent further trauma AMBU Bag 2 Device used for artificial ventilation (Sizes-adult, paediatrics) CPR Mannequin 1 Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Central Line 2 Cannula used to place in the central veins Used to open a port in blood vessel for medication and sampl withdrawal (Sizes-16,18,20,22,24)			Cylindrical Flexible rubber
CPR Board 1 Board used for to make surface firm for effective CPR Used for lifting and moving the patient in case of trauma with starps to tie the patient(orange in color with grooves cut) Spine Immobilisation Kit Splints (Rigid and Soft) AMBU Bag CPR Mannequin Defibrillator Defibrillator Central Line Cannula (diff size) Cannula (diff size) CPR Board 1 Board used for to make surface firm for effective CPR Used for lifting and moving the patient (orange in color with grooves cut) Kit used for spine immobilization with head blocks Device used for immobilizing broken bone of the body to prevent further trauma Device used for artificial ventilation (Sizes-adult, paediatrics) to learn how to give cpr (Brand-Laedral) Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Cannula used to place in the central veins Used to open a port in blood vessel for medication and sample withdrawal (Sizes-16,18,20,22,24)			instrument used for making
CPR Board 1 Board used for to make surface firm for effective CPR Used for lifting and moving the patient in case of trauma with starps to tie the patient(orange in color with grooves cut) Spine Immobilisation Kit Splints (Rigid and Soft) AMBU Bag CPR Mannequin Defibrillator Defibrillator Central Line Cannula (diff size) Cannula (diff size) CPR Board 1 Board used for to make surface firm for effective CPR Used for lifting and moving the patient (orange in color with grooves cut) Kit used for spine immobilization with head blocks Device used for immobilizing broken bone of the body to prevent further trauma Device used for artificial ventilation (Sizes-adult, paediatrics) to learn how to give cpr (Brand-Laedral) Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Cannula used to place in the central veins Used to open a port in blood vessel for medication and sample withdrawal (Sizes-16,18,20,22,24)	Bougie	2	passages
CPR Board 1 firm for effective CPR Used for lifting and moving the patient in case of trauma with starps to tie the patient(orange in color with grooves cut) Spine Immobilisation Kit 2 Splints (Rigid and Soft) 4 Device used for immobilizing broken bone of the body to prevent further trauma Device used for artificial ventilation (Sizes-adult, paediatrics) to learn how to give cpr (Brand-Laedral) Defibrillator 1 Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Central Line 2 Cannula used to place in the central veins Cannula (diff size) 5 (Sizes-16,18,20,22,24)	Dougle		Board used for to make surface
Spine Board (with straps) Spine Immobilisation Kit Splints (Rigid and Soft) AMBU Bag CPR Mannequin Defibrillator Central Line Cannula (diff size) Cannula (diff size) Cannula (diff size) Spine Board (with starps to tie the patient (orange in color with grooves cut) Kit used for spine immobilization with head blocks Kit used for spine immobilization with head blocks Color with grooves cut) Kit used for spine immobilization with head blocks Device used for immobilizing broken bone of the body to prevent further trauma Device used for artificial ventilation (Sizes-adult, paediatrics) to learn how to give cpr (Brand-Laedral) Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Cannula used to place in the central veins Used to open a port in blood vessel for medication and sample withdrawal (Sizes-16,18,20,22,24)	CPR Board	1	
Spine Board (with straps) Spine Immobilisation Kit Splints (Rigid and Soft) AMBU Bag CPR Mannequin Defibrillator Central Line Cannula (diff size) Cannula (diff size) Cannula (diff size) Cannula (diff size) Spine Immobilisation Kit used for spine immobilization with head blocks with the sad blocks at the patient (orange in color with grooves cut) Kit used for spine immobilization with head blocks Cell Device used for immobilizing broken bone of the body to prevent further trauma Device used for artificial ventilation (Sizes-adult, paediatrics) to learn how to give cpr (Brand-Laedral) Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Cannula used to place in the central veins Used to open a port in blood vessel for medication and sample withdrawal (Sizes-16,18,20,22,24)	CI R Board		Used for lifting and moving the
Spine Board (with straps) Spine Immobilisation Kit Splints (Rigid and Soft) AMBU Bag CPR Mannequin Defibrillator Defibrillator Central Line Cannula (diff size) Spine Immobilisation (Sizes-16,18,20,22,24) Spine Immobilisation with grooves cut) Kit used for spine immobilization with head blocks Device used for immobilizing broken bone of the body to prevent further trauma Device used for artificial ventilation (Sizes-adult, paediatrics) to learn how to give cpr (Brand-Laedral) Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Cannula used to place in the central veins Used to open a port in blood vessel for medication and sample withdrawal (Sizes-16,18,20,22,24)			patient in case of trauma with
Spine Immobilisation Kit Splints (Rigid and Soft) AMBU Bag CPR Mannequin Defibrillator Defibrillator Central Line Cannula (diff size) Spine Immobilisation Kit used for spine immobilization with head blocks Device used for immobilizing broken bone of the body to prevent further trauma Device used for artificial ventilation (Sizes-adult, paediatrics) to learn how to give cpr (Brand-Laedral) Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Cannula used to place in the central veins Used to open a port in blood vessel for medication and sampl withdrawal (Sizes-16,18,20,22,24)	Spine Board (with		
Spine Immobilisation Kit 2 Device used for immobilizing broken bone of the body to prevent further trauma Device used for artificial ventilation (Sizes-adult, paediatrics) to learn how to give cpr (Brand- Laedral) Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Central Line Central Line Cannula (diff size) Spine Kit used for spine immobilization with head blocks Device used for immobilizing broken bone of the body to prevent further trauma Device used for artificial ventilation (Sizes-adult, paediatrics) to learn how to give cpr (Brand- Laedral) Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Cannula used to place in the central veins Used to open a port in blood vessel for medication and sample withdrawal (Sizes-16,18,20,22,24)		1	
Immobilisation Kit 2 Device used for immobilizing broken bone of the body to prevent further trauma Device used for artificial ventilation (Sizes-adult, paediatrics) CPR Mannequin Defibrillator Defibrillator 1 Device used for artificial ventilation (Sizes-adult, paediatrics) to learn how to give cpr (Brand-Laedral) Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Cannula used to place in the central veins Used to open a port in blood vessel for medication and sample withdrawal (Sizes-16,18,20,22,24)		•	
Splints (Rigid and Soft) Splints (Rigid and Soft) AMBU Bag CPR Mannequin Defibrillator Defibrillator Central Line Cannula (diff Size) Cannula (diff Size) Splints (Rigid and Soft) Device used for immobilizing broken bone of the body to prevent further trauma Device used for artificial ventilation (Sizes-adult, paediatrics) to learn how to give cpr (Brand-Laedral) Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Cannula used to place in the central veins Used to open a port in blood vessel for medication and sample withdrawal (Sizes-16,18,20,22,24)			
Splints (Rigid and Soft) AMBU Bag CPR Mannequin Device used for immobilizing broken bone of the body to prevent further trauma Device used for artificial ventilation (Sizes-adult, paediatrics) to learn how to give cpr (Brand-Laedral) Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Central Line Central Line Central Line Cannula used to place in the central veins Used to open a port in blood vessel for medication and sample withdrawal (Sizes-16,18,20,22,24)		2	
Splints (Rigid and Soft) 4 broken bone of the body to prevent further trauma Device used for artificial ventilation (Sizes-adult, paediatrics) to learn how to give cpr (Brand-Laedral) Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Central Line Central Line Cannula used to place in the central veins Used to open a port in blood vessel for medication and sample withdrawal (Sizes-16,18,20,22,24)	Kit		Device used for immobilizing
Soft) 4 prevent further trauma Device used for artificial ventilation (Sizes-adult, paediatrics) to learn how to give cpr (Brand-Laedral) Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Central Line Central Line Cannula (diff size) Soft) 4 prevent further trauma Device used for artificial ventilation (Sizes-adult, paediatrics) to learn how to give cpr (Brand-Laedral) Caedral Cannula used to stun myocardium and deliver shock (biphasic automated defibrillator) Cannula used to place in the central veins Used to open a port in blood vessel for medication and sample withdrawal (Sizes-16,18,20,22,24)	Calinta (Digid and		broken bone of the body to
Device used for artificial ventilation (Sizes-adult, paediatrics) to learn how to give cpr (Brand-Laedral) Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Central Line Central Line Cannula (diff size) Size) Device used for artificial ventilation (Sizes-adult, paediatrics) to learn how to give cpr (Brand-Laedral) Caedral Cannula used to stun myocardium and deliver shock (biphasic automated defibrillator) Cannula used to place in the central veins Used to open a port in blood vessel for medication and sample withdrawal (Sizes-16,18,20,22,24)		4	
ventilation (Sizes-adult, paediatrics) to learn how to give cpr (Brand- Laedral) Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Cannula used to place in the central Line Central Line 2 Cannula (diff size) 5 ventilation (Sizes-adult, paediatrics) to learn how to give cpr (Brand- Laedral) Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Cannula used to place in the central veins Used to open a port in blood vessel for medication and sample withdrawal (Sizes-16,18,20,22,24)	3011)	+	
AMBU Bag CPR Mannequin Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Central Line Central Line Cannula (diff size) Cannula (diff size) Cannula used to place in the withdrawal (Sizes-16,18,20,22,24)			
to learn how to give cpr (Brand-Laedral) Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Central Line Central Line Central Line Cannula (diff size) Size To learn how to give cpr (Brand-Laedral) Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Cannula used to place in the central veins Used to open a port in blood vessel for medication and sample withdrawal (Sizes-16,18,20,22,24)	AMBII Bag	2	A Discourse of the Control of the Co
CPR Mannequin Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Cannula used to place in the central veins Used to open a port in blood vessel for medication and sample withdrawal (Sizes-16,18,20,22,24)	AMDO Dag		
Device used to stun myocardium and deliver shock (biphasic automated defibrillator) Cannula used to place in the central veins Used to open a port in blood vessel for medication and sample withdrawal (Sizes-16,18,20,22,24)	CDD Mannaquin	1	Laedral)
Defibrillator 1 automated defibrillator) Cannula used to place in the central veins Used to open a port in blood vessel for medication and sampl withdrawal (Sizes-16,18,20,22,24)	CFR Mainiequii	+	
Defibrillator 1 automated defibrillator) Cannula used to place in the central veins Used to open a port in blood vessel for medication and sample withdrawal (Sizes-16,18,20,22,24) Size 1 (Sizes-16,18,20,22,24)			and deliver shock (biphasic
Central Line 2 Cannula used to place in the central veins Used to open a port in blood vessel for medication and sample withdrawal (Sizes-16,18,20,22,24) 5 Cannula used to place in the central veins Used to open a port in blood vessel for medication and sample withdrawal	Defibrillator	1	
Central Line 2 central veins Used to open a port in blood vessel for medication and sample withdrawal (Sizes-16,18,20,22,24) 5 (Sizes-16,18,20,22,24)	Denomiator	 	
Cannula (diff size) Used to open a port in blood vessel for medication and sample withdrawal (Sizes-16,18,20,22,24)	Central Line	2	
Cannula (diff size) vessel for medication and sample withdrawal (Sizes-16,18,20,22,24)	Cellulai Lille		
Cannula (diff size) 5 (Sizes-16,18,20,22,24)			vessel for medication and sample
size) 5 (Sizes-16,18,20,22,24)	Cannula (diff		
6: 5 110 120 150 1		5	
	Syringes	4	Sizes-5ml,10ml,20ml,50ml

		To help control bleeding, and
		immobilizing, dressing
Bandages	4	Sizes-4inch,6inch
Foley		Tube placed in the ureter for
Cathether(diff		passing the urine
	4	Size-14,16,18,20
size)		Drugs used in emergency
,		situation(Atropine,NORAD,Epin
		ephrine,
		Hydrocortisone,
		Pantop, buscopan,
		Drotin, labetolol, midazolam, calci
		um gluconate, Trenaxa, Normal
		saline etc)
	16	January 1117
Emergency Drugs	15	Device used for ventilation
		where ET tube placement is not
	1	feasible
Combitube		Device used for ventilation
		Device used for ventuation
Laryngeal Mask		where ET tube placement is not
Airway	11	feasible C Spine
		To immobilize the C-Spine
Cervical Collar	2	(Soft&Rigid)
		Instrument to put in buccal cavity
Oropharyngeal		for ventilation and preventing
Airway	5	toungue bite (all size)
7 th may		To remove excess fluid mucus
Suction Catheter	2	from the body
Suction Cameter		Basic emergency kit with all the
		medical
Partie Language		supply(bandaid,bandages,microp
		ore,
		syringes
	1	diameter de la constant de la consta
First Aid Kit	- 1	



School of Nursing Sciences and Allied Health Department of Paramedical Sciences Jamia Hamdard

List of Equipment of Dialysis under Rs.25,000/-

S.no.	Name of Equipment	Quantity	Specifications
1.	Arrow Dialysis single lumen catheter	01	Single lumen catheter along with guide wire, dilator, introducer needle.
2.	Arrow Dialysis double lumen catheter	01	Double lumen catheter along with guide wire, dilator, introducer needle
3.	Fistula Needle 16 gauze (02) 17 gauze (04) 18 gauze (02)	08	Fistula needle with backeye of different sizes.
4.	Normal Saline 500ml 250ml	10	For demo purpose.
5.	I.V Sets	10	For demo purpose
6.	Betadine	02	For demo purpose

7.	Sprit	02	For demo purpose.
8.	Sterile gauze packets	20	For demo purpose.
9.	Syringes and needles (2ml, 5ml, 10ml)	10	For demo purpose
10.	Sterile HD Tray	05	For demo
			purpose.